

Maintenance of *Arabidopsis thaliana* T87 cell suspension culture

Introduction

T87 cell suspension culture was obtained from a seedling of *Arabidopsis thaliana* L. (Heynh.) ecotype Columbia (Axelos *et al.* 1992). T87 cell suspension culture is composed of small, near-uniform clumps of cells and has a green color (Figure 1). T87 cells have been maintained in Jouanneau and Péaud-Lenoël (JPL) medium under continuous light at 22°C and subcultured at two-week intervals.

Materials

I. Stock Solutions and Chemicals

A) JPL A' (1 L)		
KNO ₃		65.5 g
CaCl ₂ ·2H ₂ O		4.4 g
MgSO ₄ ·7H ₂ O		3.7 g
KH ₂ PO ₄		1.7 g
B) JPL B (1 L)		
H ₃ BO ₃		6.2 g
MnSO ₄ ·4H ₂ O		22.3 g
ZnSO ₄ ·7H ₂ O		10.6 g
KI		0.83 g
Na ₂ MoO ₄ ·2H ₂ O		0.25 g
CoCl ₂ ·6H ₂ O		0.025 g
CuSO ₄ ·5H ₂ O		0.025 g
C) JPL C (1 L)		
FeSO ₄ ·7H ₂ O		2,780 mg
Na ₂ -EDTA		3,730 mg
D) JPL D (1 L)		
<i>myo</i> -Inositol		10 g
Glycine		0.2 g
E) JPL VT (100 mL)		
Nicotinic acid		50 mg
Pyridoxine·HCl		50 mg
Thiamine·HCl		40 mg
F) JPL P (100 mL)		
200 mM KH ₂ PO ₄		19.5 ml
200 mM Na ₂ HPO ₄		30.5 ml
H ₂ O		50 ml

- G) 1 mM NAA (100 mL)
 1-Naphtaleneacetic acid 18.62 mg
- H) Casein hydrolysate, vitamin-free
 Casamino acids vitamin assay, Difco (228820)

II. Glassware and Stainless Sieves

- A) Erlenmeyer flask (300 ml), capped with two layers of aluminum foil
 - B) Pipette, large tip opening (10 ml), and a bulb
 - C) Stainless sieve (diameter, 5 cm; pore size, 1 mm) set on a tall beaker (200 ml), capped with two layers of aluminum foil
- All are sterilized by autoclaving at 121°C for 20 min.

III. Preparation of JPL Medium

- A) Prepare three solutions as follows:

- A)-1 JPL mineral solution (1 L)

JPL A'	37.5 ml
JPL B	0.375 ml
JPL C	2.5 ml

 Adjust pH to 5.7 with 0.2 N KOH.

- A)-2 JPL organic solution (100 mL)

Casein hydrolysate	0.1 g
JPL D	10 ml
JPL VT	1 ml

 Adjust pH to 5.7 with 0.2 N HCl.

- A)-3 JPL sucrose solution (100 mL)

Sucrose	15 g
JPL P	1 ml
1 mM NAA	1 ml

- B) Autoclave these solutions at 121°C for 20 min.
- C) Add 800 ml of JPL mineral solution, 100 ml of JPL organic solution, and 100 ml of JPL sucrose solution aseptically.
- D) Pour 80 ml of JPL medium into a sterile flask.

Methods

- I. Filter a two-week-old cell suspension through a stainless sieve (Figure 2).
- II. Agitate the filtrate well and transfer 2 ml of cell suspension to 80 ml of fresh JPL medium with a pipette.
- III. Incubate cell cultures on a rotary shaker at 120 rpm under continuous light (40–100 $\mu\text{E/s/m}^2$) at 22°C.

Notes

- I. T87 cells may be maintained in other media. We have cultured T87 cells in JPL medium as described in the original report by Axelos *et al.* (1992). Yamada *et al.* (2004) reported that T87 cells were well grown in Gamborg's B5 medium. When T87 cells are cultured in the medium other than JPL medium, it is necessary to optimize the subculture method, such as subculture interval, amount of cells transferred, and the pore size of the sieve.
- II. T87 cells should be transferred to fresh medium immediately after arrival. We send T87 cells on semi-solid JPL medium in 250-ml disposable flasks. T87 cells are collected from the medium by a spatula and transferred to Erlenmeyer flasks containing fresh liquid medium. The cell suspension culture should be established from small scale (e.g. 20 ml of medium in a 100-ml flask), because T87 cells may have been damaged during transport.
- III. The amount of living cells subcultured to fresh medium is important for stable maintenance of cell suspension cultures. Higher-density culture causes overgrowth after the usual culture period, while lower-density culture suppresses the cell division. Occasionally most of T87-cell clumps develop into large aggregates, which cause a decrease in the number of cells passed through a 1-mm sieve. In such cases, it should be confirmed that an adequate amount of cells are transferred to each fresh medium during subculturing.
- IV. In our culture condition, T87 cells increased around 50-fold on the day 14 of culture (Figure 3). In addition, T87 cells also proliferated in the dark.

References

- Axelos M, Curic C, Mazzolini L, Bardet C, Lescure B (1992) A protocol for transient gene expression in *Arabidopsis thaliana* protoplasts isolated from cell suspension cultures. *Plant Physiology and Biochemistry* 30: 123-128
- Yamada H, Koizumi N, Nakamichi N, Kiba T, Yamashino T, Mizuno T (2004) Rapid response of *Arabidopsis* T87 cultured cells to cytokinin through His-to-Asp phosphorelay signal transduction. *Bioscience, Biotechnology, and Biochemistry* 68: 1966-1976





Figure 1 *Arabidopsis thaliana* T87 cell suspension culture

- (A) Two-week-old cell suspension culture.
- (B) Cell clumps. The size of cell clumps are varied, but most of them are below 1 mm.
- (C) Microscopic observation of T87 cells. The cells contain many chloroplasts.



2-week-old T87 suspension culture in a 300-ml flask.



Pass the cell suspension through a stainless sieve.

Stainless sieve (pore size, 1 mm)



Transfer 2 ml of filtrate to 80 ml of fresh JPL medium.

Figure 2 Procedure for subculturing T87 cells




Figure 3 Growth profiles of T87 cells